ATS/ERS Pulmonary Function Acceptability & Repeatability Guidelines

Forced Vital Capacity



ACCEPTABILITY

- 1. Back extrapolated volume must be <5% of the FVC or 0.150 L, whichever is greater.
- 2. The volume/time curve shows no change in volume for ≥1 second (subject has a plateau).
- For subjects 7-9 years old, the subject has tried to exhale for ≥3 seconds OR for subjects ≥10 years old, the subject has tried to exhale for ≥6 seconds (>15 seconds rarely changes clinical decisions).
- 4. There should be no cough detected in the first second of exhalation which could affect FEV1.
- 5. There should be no leaks, Valsalva maneuver (glottis closure) or obstruction of mouthpiece
- 6. An adequate test requires a minimum of three acceptable FVC maneuvers with two of them meeting repeatability criteria.

ACCEPTABILITY PRESCHOOL CHILDREN (<6 YEARS OLD)

- 1. Back extrapolated volume of <80 mL or 12.5%
- 2. Cessation of flow occurs at 10% or less of Peak Flow.
- 3. Expiratory time should be reported but not used to exclude maneuvers.
- 4. At least two acceptable maneuvers.

Diffusing Capacity



ACCEPTABILITY

- 1. Inspired volume should be \geq 90% of largest Vital Capacity*.
- 2. 85% of test gas inhaled in <4 seconds.
- 3. Breath Hold Time should be between 8-12 seconds.
- 4. Sample collection should be completed within 4 seconds of the start of exhalation. For RGA systems, virtual sample collection should be initiated after dead-space washout is complete.
- 5. No evidence of leaks, or Valsalva or Mueller maneuvers during lockout.
- 6. At least 4 minutes between tests to allow an adequate elimination of test gas from the lungs for classical systems. For RGA systems, tracer gas level at end-exhalation must be ≤2% of the tracer gas concentration in the test gas.

REPEATABILITY

- 1. The difference between the largest and next largest FVC is ≤ 0.150 L.
- 2. The difference between the largest and next largest FEV1 is ≤ 0.150 L.
- 3. If FVC is ≤ 1.0 L, these values are reduced from 0.150 L to 0.100 L.
- 5. The largest values for FVC, FEV1 and Peak Flow, from three acceptable maneuvers are reported.

REPEATABILITY PRESCHOOL CHILDREN (<6 YEARS OLD)

- 1. The difference between the largest and next largest FVC is ≤ 0.1 L.
- 2. The difference between the largest and next largest FEV1 is \leq 0.1 L.

BRONCHODILATOR RESPONSE

- 1. Efforts can be recorded ≥10 minutes and up to 15 minutes after giving short-acting B2-agonists, and 30 minutes after giving short-acting anticholinergic agents.
- 2. A >12% and 200 mL change in FEV1 compared to baseline suggests a "significant" bronchodilatation OR a >12% and 200 mL change in FVC compared to baseline suggests a "significant" bronchodilatation.

REPEATABILITY

- 1. DLco values should be within 2 mL/min/mmHg (0.67 mmol/min/kPa) of each other.
- 2. The average of at least two acceptable tests that meet repeatability should be reported.

NOTE: Adjustments of DLco for Hb, COHB and altitude should be considered.

*A maneuver with an inspired volume of \geq 85% of largest vital capacity may be deemed acceptable if the VA is within 200mL or 5% (whichever is greater) of the largest VA from other acceptable maneuvers.

TGV/FRCpleth



ACCEPTABILITY

- 1. Closed shutter panting frequency between 0.5 and 1.0 Hz (30-60/min).
- 2. Patient's cheeks are to be supported by both hands and the subject should breathe quietly until a stable end-expiratory level is achieved (usually 3-10 tidal breaths) before closing the shutter.
- 3. A series of 3-5 technically satisfactory panting maneuvers should be recorded.

NOTE: It is recommended to perform an SVC immediately after the shutter reopens in order to accurately calculate Total Lung Capacity.

Nitrogen Washout



ACCEPTABILITY

- 1. N2 concentration should be <1.5% for at least three successive breaths before ending test.
- 2. A change in inspired N2 of >1% or sudden large increases in expiratory N2 concentrations may indicate a leak.
- 3. At least one technically satisfactory measurement should be obtained.

REPEATABILITY

- 1. At least three TGV (FRCpleth) values that agree within 5% (the difference between the highest and lowest value divided by the mean).
- 2. The average value should be reported.

REPEATABILITY

1. If more than one measurement of FRCN₂ is made, the value reported should be the average of technically acceptable results that agree within 10%.

NOTE: If more than one washout is performed, a waiting period of \geq 15 minutes is recommended between trials. Patients with severe COPD should wait \geq 1 hour between trials.

For additional information on pulmonary diagnostics, please consult the ATS/ERS guidelines:

www.thoracic.org || www.ers-education.org

references

- ATS/ERS Task Force: Standardisation of the measurement of lung volumes, Eur Respir J 2005; 26: 511-522
- ATS/ERS Task Force: Standardisation of spirometry, Eur Respir J 2005; 26: 319-338 ATS/ERS Standards for Single-Breath Carbon Monoxide Uptake in the Lung, Eur Respir J 2017; 49; 1600016
- ATS/ERS Statement: Pulmonary Function Testing in Preschool Children, AM J Respir Crit Care Med; Vol 175: 1304-1345, 2007

MGC DIAGNOSTICS CORPORATION through its subsidiary Medical Graphics Corporation

350 Oak Grove Parkway St. Paul, Minnesota USA 55127-8599

www.mgcdiagnostics.com

© 2017 MGC Diagnostics Corporation or one of its affiliates.

All rights reserved.

