RESMON PRO

V3

TESTING AND RESULTS EVALUATION



WHAT IS THE RESMON PRO FULL?

It is a device based on the Forced Oscillation Technique (FOT), offering a complete functional assessment of the respiratory system, through simple measurements performed at tidal breathing.



MGC Diagnostics Corporation Manufactured by ResTech distributed exclusively by Medical Graphics Corporation.

350 Oak Grove Parkway St. Paul. Minnesota, USA 55127-8599 mgcdiagnostics.com

© 2020 MGC Diagnostics Corporation or one of its affiliates. All rights reserved

Designed, developed, manufactured by:



FDA CLEARED

RESTECH is an ISO 13485, ISO 9001 and MDSAP certified company.

CE

0051

Part# 035008-001 RevA

PROPER POSTURE DURING THE TEST

Maintain a sitting position, with a straight back

while leaning against the seat back and with a

mouthpiece, keeping the tongue below,

3. Support from behind the patient's cheeks

and the soft tissue under chin during the

sides to obtain a relaxed shoulders posture,

test with the patient arms falling on the

see figure below (suggested technique).

Alternatively, the patient may support his/ her cheeks ensuring that the elbows are

slightly detached from the chest.

SUGGESTED TESTING TECHNIQUE

ALTERNATIVE TECHNIQUE

IMPORTANT NOTE

Validity of results depends from good data quality and correct

breath-reject algorithms will minimize artifacts such as glottis

closure, coughs, and irregular breathing.

testing procedure, patient head/neck, relaxed shoulders position,

and supported cheeks. The Resmon Pro's sophisticated software

relaxed, slightly overextended neck.


2. Breathe in a relaxed way through the

1. Wear the noseclip.


avoiding leaks.

TEST EVALUATION

1. PRESENCE OF RESPIRATORY IMPAIRMENT AND REVERSIBILITY




Resistance (Rrs) graphs for Inspiratory, Expiratory and Total inspiratory cycle parameters at the lowest measured frequencies (for adult and pediatrics). Predicted dotted line and ULN (Upper Limit of Normality).





RRS > ULN and/or XRS < LLN are indicative of an anomaly in respiratory mechanics.

2. LOCALIZATION



and Reactance (Xrs) do not present anomalies (Rrs < ULN and Xrs > LLN).






Resistance (Rrs) is above its upper limit of normality (Rrs > ULN) and Reactance (Xrs) does not present anomalies



('T)]

z[CI

3. TIDAL EXPIRATORY FLOW LIMITATION, **ΔXRS** INDEX



Both Resistance (Rrs)


NORMAL

CENTRAL OBSTRUCTION (Xrs < LLN) for diseases affecting central airways.



Reactance (Xrs) graphs for Inspiratory, Expiratory and Total inspiratory cycle parameters at the lowest measured frequencies (for adult and pediatrics). Predicted dotted line and LLN (Lower Limit of Normality)

Differences between tests that are above those expected in a reference healthy population are highlighted in red.



5, 11, 19 Hz

PERIPHERAL DISEASE

Resistance (Rrs) does not present anomalies (Rrs < ULN), Reactance (Xrs) is below its lower limit of normality (Xrs < LLN), for possible small airway obstruction, excluded alveoli, disomogeneity of ventilation, or possible restriction.

SEVERE OBSTRUCTIVE DISEASE

Both Resistance (Rrs) is above its upper limit of normality (Rrs > ULN) and Reactance (Xrs) is below its lower limit of normality (Xrs < LLN). Resistance (Rrs) tends to decrease at higher frequencies (i.e. severe asthma, severe COPD).

5, 11, 19 Hz) (5 Hz

∆Xrs is the patented index of expiratory flow limitation during tidal breathing.*

∆Xrs> 2.8 → LIMITATION

* Dellacà et al Eur Resp. 12004 Fur Respir J 2007