



INTERPRETING FeNO READINGS

ATS/ERS CLINICAL GUIDELINES SUMMARY FOR INTERPRETING FeNO LEVELS


Diagnosis using the Fenom Flo Device			
FeNO (ppb) Levels	LOW <25ppb (<20ppb in children)	INTERMEDIATE 25-30ppb (20-35 ppb in children)	HIGH >50ppb (>35ppb in children) or rise in FeNO of >40% from previously stable levels
Symptomatic (chronic cough, and/or wheeze and/or shortness of breath during past 6 weeks)	Eosinophilic airway inflammation unlikely Alternative diagnosis Unlikely to benefit from ICS	Be cautious Evaluate clinical context Monitor change in FeNO over time	Eosinophilic airway inflammation present Likely to benefit from ICS

Alternative considerations (if Allergic Asthma has been dismissed)			
Non-Allergic Asthma	Chronic Cough	Vocal Cord Disfunction	Gerd

Monitoring (in patients with diagnosed astham) using the Fenom Flo device			
FeNO (ppb) Levels	LOW <25ppb (<20ppb in children)	INTERMEDIATE 25-30ppb (20-35 ppb in children)	HIGH >50ppb (>35ppb in children) or rise in FeNO of >40% from previously stable levels
Symptomatic (chronic cough and/or wheeze)	Possible alternative diagnosis Unlikely to benefit from increase in ICS	Persistent allergen exposure Inadequate ICS dose Poor adherence Sterioid resistance	Persistent allergen exposure Poor adherence or inhaler technique Inadequate ICS dose Risk of Exacerbation
Symptoms Absent	Adequate ICS dose Good adherence ICS taper	Adequate ICS dosing Good adherence Monitor change in FeNO	Steroid resistance ICS withdrawal or dose reduction may result in relapse

REFERENCES:

1. Dweik RA et al. An official ATS clinical practice guideline: interpretation of exhaled nitric oxide levels (FeNO) for clinical applications. Am J Respir Crit Care Med. 2011;184(5):602-15.
2. Saito J, Gibeon D, Macedo P, et al. Domiciliary diurnal variation of exhaled nitric oxide fraction for asthma control. The European Respiratory Journal. 2014 Feb;43(2):474-484. DOI: 10.1183/09031936.00048513. PMID: 23949962.



MGC DIAGNOSTICS CORPORATION, through its subsidiary Medisoft S.A.
350 Oak Grove Parkway St. Paul, Minnesota USA 55127-8599

 Medisoft S.A. P.A.E de Sorinnes, Rue du Clairon 5
5503 Sorinnes BELGIUM

© 2024 MGC Diagnostics Corporation or one of its affiliates. All rights reserved.

All specifications subject to change without notice. Products may vary from those illustrated.

MGC Diagnostics and its affiliates are equal opportunity/affirmative action employers committed to cultural diversity in the workforce.

Part# 060182-001 RevA

