

# BreezeSuite™ Cardiorespiratory Diagnostic Software

Version 8.6  
Installers Guide



## Contact

### **MGC Diagnostics Corporation through its subsidiary Medical Graphics Corporation**

350 Oak Grove Parkway  
St. Paul, Minnesota USA 55127-8599  
T +1 651.484.4874  
800.950.5597  
F +651.484.8941



# CONTENTS

Introduction.....	3
Purpose of this Manual.....	3
Professional Installation Service .....	3
Standalone Installation Steps .....	4
Installation .....	4
BreezeSuite MultiUser Software Installation Steps.....	4
Installation .....	4
Standalone Workstation Upgrade Steps .....	5
Upgrade .....	5
BreezeSuite MultiUser Software Upgrade Steps .....	5
Upgrade .....	5

# INTRODUCTION

## Purpose of this Manual

To provide guidance for successful installation or upgrade of:

- BreezeSuite cardiorespiratory diagnostic software and BreezeReview™ physician review software (standalone)
- BreezeSuite cardiorespiratory diagnostic software and BreezeSuite MultiUser™ networking software

## Professional Installation Service

MGCD Diagnostics offers professional installation services (service fees apply).

Professional service installation benefits include:

- Remote or onsite installation service
- Scheduled and reliable software installation for minimal downtime
- Software configuration assistance
- Report customization assistance
- Predicted set configuration assistance
- Diagnostics testing screen customization assistance
- Application security configuration guidance
- Data conversion assistance
- Device communication verification


Please contact MGC Diagnostics to discuss the installation service options available.

Phone: 800-333-4137

Email: [Support@MGCDiagnostics.com](mailto:Support@MGCDiagnostics.com)

## STANDALONE INSTALLATION STEPS

Follow the installation process details on this page if you are installing software on a PC without connecting to a shared/centralized database.

 **Important:** MGC Diagnostics strongly recommends investing in remote or onsite installation assistance where a qualified and experienced Professional Service Representative will ensure the installation is completed successfully with minimal impact on lab workflow.

### Installation


For successful standalone software installation, perform the following steps:


- Print a copy of this page to follow during the installation.
- Follow the instructions detailed in the manuals and guides listed below in the order provided.
- Check off each manual or guide listed after each task has been completed.

- ☐ 1. Standalone Workstation Installation Manual
- ☐ 2. Security Administrators Guide
- ☐ 3. Workstation Configuration Guide

## BREEZESUITE MULTIUSER SOFTWARE INSTALLATION STEPS

Follow the installation process details on this page if you install software on a computer and connect to a Network SQL Server Database.

 **Note:** Software installation requires access to a server with SQL Server installed and network connectivity between the Computer(s) and SQL Server. Software installation requires investment in the BreezeSuite MultiUser networking software option.

 **Important:** MGC Diagnostics strongly recommends investing in remote or onsite installation assistance where a qualified and experienced Professional Service Representative will ensure the installation is completed successfully with minimal impact on lab workflow.

### Installation


For successful standalone software installation, perform the following steps:

- Print a copy of this page to follow during the installation.
- Follow the instructions detailed in the manuals and guides listed below in the order provided.
- Check off each manual or guide listed after each task has been completed.

- ☐ 1. MultiUser Server Installation Manual
- ☐ 2. MultiUser Workstation Installation Manual
- ☐ 3. Security Administrators Guide
- ☐ 4. Workstation Configuration Guide

## STANDALONE WORKSTATION UPGRADE STEPS

Follow the process details on this page if you are upgrading from BreezeSuite 8.1 or higher software on a client workstation without connecting to a Network SQL Server Database.

 **Important:** MGC Diagnostics strongly recommends investing in remote or onsite installation assistance where a qualified and experienced Professional Service Representative will ensure the installation is completed successfully with minimal impact on lab workflow.

### Upgrade

For successful standalone software installation, perform the following steps:


- Print a copy of this page to follow during the installation.
- Follow the instructions detailed in the manuals and guides listed below in the order provided.
- Check off each manual or guide listed after each task has been completed.

- ☐ 1. Standalone Workstation Upgrade Manual
- ☐ 2. Security Administrators Guide

## BREEZESUITE MULTIUSER SOFTWARE UPGRADE STEPS

Follow the process details on this page if you are upgrading from BreezeSuite 8.1 or higher software on a client computer connected to a Network SQL Server Database

 **Note:** The software upgrade requires an upgrade of the databases on the Network SQL Server

 **Important:** MGC Diagnostics strongly recommends having a support agreement in place prior to installing software. A support agreement is required for installation assistance and ongoing support of software.

### Upgrade

For successful standalone software installation, perform the following steps:

- Print a copy of this page to follow during the installation.
- Follow the instructions detailed in the manuals and guides listed below in the order provided.
- Check off each manual or guide listed after each task has been completed.

- ☐ 1. MultiUser Server Upgrade Manual
- ☐ 2. MultiUser Workstation Upgrade Manual
- ☐ 3. Security Administrators Guide