



RAM Treadmills


Stress Testing Treadmills



Clinical 870A



Sport 870S



Athlete 870C

CHOOSE YOUR MODEL

[Model: Clinical 870A]

Footprint: 31.5"W x 80.7"L x 44.1"H
(80 x 205 x 112 cm)
Height is with treadmill at zero inclination

Running Surface: 19.7"W x 59"L
(50 x 150 cm)

Incline Range: 0-25%

Speed Range: 0-12 mph (0-20 km/h)

PC Console with 15.6" touchscreen

Lanyard Auto-Stop

[Model: Sport 870S]

Footprint: 35"W x 90.6"L x 45.3"H
(89 x 230 x 115 cm)
Height is with treadmill at zero inclination

Running Surface: 22.8" x 66.9"
(58 x 170 cm)

Incline Range: 0-22%

Speed Range: 0-16 mph (0-25 km/h)

PC Console with 15.6" touchscreen

Lanyard Auto-Stop

[Athlete 870C]

Footprint: 37.8"W x 99.2"L x 45.3"H
(96 x 230 x 115 cm)
Height is with treadmill at zero inclination

Running Surface: 25.6" x 74.8"
(65 x 190 cm)

Incline Range: 0-22%

Speed Range: 0-20 mph (0-32 km/h)

PC Console with 15.6" touchscreen

Lanyard Auto-Stop

PERFECT FOR CLINICAL, REHABILITATION, AND SPORTS MEDICINE

RAM treadmills are designed to meet the rigorous demands of healthcare and sports medicine settings. Our treadmills offer a range of features to ensure optimal performance and patient safety:

SUPERIOR COMFORT: Our triple-thickness running belt is engineered to provide maximum energy absorption, making every workout as comfortable as possible.

ENHANCED SAFETY: Non-slip rubber coating and parallel bars offer superior stability and grip, minimizing the risk of accidents.

UNMATCHED DURABILITY: Reinforced rigidity ensures longevity and reliable performance, even in high-traffic environments.




WHY CHOOSE RAM TREADMILLS?

- **Unmatched Italian Craftsmanship:** RAM treadmills are built with the finest materials and meticulous attention to detail, ensuring exceptional durability and a smooth running experience.
- **Advanced Technology:** RAM treadmills incorporate innovative features like variable speed control, incline adjustment, and built-in workout programs to help you achieve your fitness goals.
- **Exceptional Performance:** RAM treadmills deliver a powerful motor and a spacious running surface to accommodate users of all fitness levels.
- **Unwavering Reliability:** RAM treadmills are built to last, providing years of dependable performance in clinical or home settings.




PLUG OPTIONS:




MANUFACTURER'S WARRANTY:

1 year parts & labor, 2 years parts

[ADDITIONAL OPTIONS FOR ALL MODELS]



Structure support for safety harness with auto-stop



Support chair movable and removeable



Pediatric Hand Rail

ALSO AVAILABLE:
Belt for HR Monitoring with Bluetooth® technology



Lanyard auto-stop

MGC DIAGNOSTICS CORPORATION, manufactured by Medisoft RAM Italia Srl, distributed by Medical Graphics Corporation
350 Oak Grove Parkway St. Paul, Minnesota USA 55127-8599

© 2024 MGC Diagnostics Corporation or one of its affiliates. All rights reserved.
All specifications subject to change without notice. Products may vary from those illustrated.
MGC Diagnostics and its affiliates are equal opportunity/affirmative action employers committed to cultural diversity in the workforce.
Part# 061384-001 RevB