



MGC DIAGNOSTICS®

through its subsidiary medisoft®

FeNO+ PFS

Configurator



www.mgcdiagnostics.com

FeNO+ PFS



● FeNO+ multi-flow (100-240 V, 50-60 Hz) (bronchial and alveolar exhaled NO) Ref: FE-00-01

Important note A FeNO+ always includes the ExpAir software with:

- Exhaled Nitric oxide (bronchial)
- Exhaled Nitric oxide (alveolar)

The standard module is vertically mounted.

Start with the base module and follow this path :

- 1 Options
- 2 Device mounting
- 3 Computer
- 4 Trolley

1 Options

○ Nasal NO measurement Ref: FE-00-02

○ Spirometry Ref: FE-00-04

- Includes :
- Resting spirometry (VC)
 - Forced spirometry (FVC)
 - Maximum voluntary ventilation (MW)

- Minute tidal ventilation (Vmin)
- Reversibility (pre/post)
- Bronchial challenge, multi-level spirometry

2 Device mounting

● Left vertical mounting on trolley Ref: FE-00-05

● Right vertical mounting on trolley Ref: FE-00-06

3 Computer See IT Hardware configurator

4 Trolley See Trolley configurator

All the options listed are available depending on the products and may be unavailable within some countries.