



**MGC  
DIAGNOSTICS®**

### Platinum Elite® body plethysmograph

diameter: 41.5 in 105.4 cm

circumference: 113 in 287 cm

height: 69.5 in 176.5 cm

system weight: 350 lbs 158 kgs

internal volume: 726 L

- Safe, simple “turn and sit” entry with stable, no-twist seat
- Zero-clearance door provides better use of available lab space
- Largest usable patient space with smallest interior volume  
(no functional patient weight limit)
- Fits through 32 in ( 81.3 cm) door