



# RAM Treadmills

Stress Testing Treadmills



Model: Clinical 870A


RAM treadmills are perfectly adapted for use in clinical, rehabilitation and sports medicine applications.

## FEATURES

- Triple thickness running belt is designed to guarantee maximum energy absorption making use as comfortable as possible.
- Non-slip rubber coating and parallel bars for maximum patient safety.
- Reinforced rigidity of the unit results in increased mechanical strength.
- 3 treadmill models

## WHY CHOOSE RAM TREADMILLS?

- **Unmatched Italian Craftsmanship:** RAM treadmills are built with the finest materials and meticulous attention to detail, ensuring exceptional durability and a smooth running experience.
- **Advanced Technology:** RAM treadmills incorporate innovative features like variable speed control, incline adjustment, and built-in workout programs to help you achieve your fitness goals.
- **Exceptional Performance:** RAM treadmills deliver a powerful motor and a spacious running surface to accommodate users of all fitness levels.
- **Unwavering Reliability:** RAM treadmills are built to last, providing years of dependable performance in clinical or home settings.



## SPECIFICATIONS

Model: Clinical 870A

### [ Standard ]

19.7" x 59" (50 x 150 cm) running surface
Incline Range: 0-25%
Speed Range: 0-12 mph (0-20 km/h)
PC Console with 15.6" touchscreen
Lanyard Auto-Stop

Model: Sport 870S

### [ Standard ]

22.8" x 66.9" (58 x 170 cm) running surface
Incline Range: 0-22%
Speed Range: 0-16 mph (0-25 km/h)
PC Console with 15.6" touchscreen
Lanyard Auto-Stop

Model: Athlete 870C


### [ Standard ]

25.6" x 74.8" (65 x 190 cm) running surface
Incline Range: 0-22%
Speed Range: 0-20 mph (0-32 km/h)
PC Console with 15.6" touchscreen
Lanyard Auto-Stop



Lanyard auto-stop

### [ Additional Options for All Models ]



Structure support for safety harness with auto-stop



Support chair movable and removeable



Pediatric Hand Rail

ALSO AVAILABLE:  
Belt for HR Monitoring with Bluetooth® technology

## MANUFACTURER'S WARRANTY:

1 year parts & labor, 2 years parts



MGC DIAGNOSTICS CORPORATION, manufactured by Medisoft RAM Italia Srl, distributed by Medical Graphics Corporation  
350 Oak Grove Parkway St. Paul, Minnesota USA 55127-8599

© 2024 MGC Diagnostics Corporation or one of its affiliates. All rights reserved.

All specifications subject to change without notice. Products may vary from those illustrated.

MGC Diagnostics and its affiliates are equal opportunity/affirmative action employers committed to cultural diversity in the workforce.

Part# 061384-001 RevA