

# MDI Holding Chambers & Spacers

MicroChamber® and MicroSpacer®



The MicroChamber® and the MicroSpacer® are devices for use with metered-dose inhaler (MDI) products. Both devices deliver aerosol medication from your MDI directly to the lungs, where it is needed the most.



Deeper Penetration,  
Faster Delivery To  
Your Lungs



Clinically Proven More  
Effective Than Larger  
Devices



Compatible With Most  
Metered Dose Inhaler  
(MDI) Products



Smaller Device For  
Easy Travel &  
Portability



Perfect For All Ages,  
Kids 6-12, Teens,  
Adults & Elderly



Quality Products,  
MADE IN THE USA!

## MicroChamber®

MicroChamber® is an aerosol **holding chamber** for use with metered-dose inhaler (MDI) products. MicroChamber® can be attached to a metered-dose inhaler (MDI) in seconds, and can be easily detached for cleaning and drying. MicroChamber® provides a chamber device that is small, convenient, effective, and less expensive than other chamber devices.



### FEATURES


- Optimizes MDI Medication benefits by delivering more medication to the lungs.
- Unique design results in faster and deeper medication delivery to the lungs.
- Clinically proven to be as effective as larger devices.
- Compatible with most mask products (22mm connector).
- Built-in air vents increase tidal volume for deeper penetration to lungs.
- Compatible with most MDI Products.
- Attaches and detaches from inhaler easily for cleaning.
- Small and portable
- Great for kids, elderly, and those needing extra time.
- Dishwasher safe, latex-free, comfortable to use.

The MicroChamber® device is made of a durable polypropylene plastic material.



# MicroSpacer®

MicroSpacer® is an aerosol **spacer** device for use with metered-dose inhaler (MDI) products. MicroSpacer® represents a technological breakthrough regarding proper airway disposition of asthma medication. MicroSpacer® reduces “impactor throat” deposition over 70%. The spacer extends the mouthpiece of a metered-dose inhaler only one inch, and weighs less than one ounce. The MicroSpacer® can be attached to a metered-dose inhaler in seconds, and can be easily detached for cleaning and drying. Each MicroSpacer® should last from four to six months before requiring replacement. MicroSpacer® outperforms the larger chamber devices regarding percentage of respirable dose delivered.



## FEATURES

- Optimizes MDI Medication benefits by delivering more medication to the lungs.
- Small size offers great portability for people on the go.
- Unique design results in faster and deeper medication delivery to the lungs.
- Clinically proven to be as effective as the larger spacer/chamber devices, but at a much less cost.
- Compatible with most MDI products.
- Attaches to MDI in seconds, and easily detaches for cleaning.
- Pocket-size design. Weighs less than 1oz and extends MDI mouthpiece by just 1 inch.
- Smaller device size ensures adequate coordination time with spraying of the inhaler.
- Ideal for children 6-12 years old, teenagers and adults on the go.
- Easy to clean, dishwasher safe, and latex free and anti-static delivery.
- Replace every 4-6 months. 100% recyclable.
- For single patient use.

The MicroChamber® device is made of a durable polypropylene plastic material.

## PART INFORMATION

MicroChamber® | Part Number: 536108-001 | Case of 50

MicroSpacer® | Part Number: 536109-001 | Case of 144

## ORDERING INFORMATION

To order, call 1.800.655.1133.

DESIGNED, DEVELOPED AND  
MANUFACTURED BY:



MGC DIAGNOSTICS CORPORATION, Manufactured by Respiratory Delivery Systems, distributed by Medical Graphics Corporation  
350 Oak Grove Parkway St. Paul, Minnesota USA 55127-8599

© 2024 MGC Diagnostics Corporation or one of its affiliates. All rights reserved.

All specifications subject to change without notice. Products may vary from those illustrated.

MGC Diagnostics and its affiliates are equal opportunity/affirmative action employers committed to cultural diversity in the workforce.

Part# 060183-001 RevA